

OEM · ODM · Medical device · Beauty device · Skin care device



FOR BEAUTIFUL SKIN FOR YOUR DREAM

HR MEDITECH
PRODUCTS CATALOG



www.hrmeditech.co.kr

Korean skin beauty equipment Leader in beauty expertise



PASSION

We aim to develop into a global leader in the beauty care market.



VISION

We lead the k-beauty culture through innovative technology development through continuous R&D investment.



CONFIDENCE

HR MEDITECH prioritizes trust with customers.

HR MEDITECH

As a global leader in the beauty & healthcare market, all employees are doing their best in the high-tech medical device and beauty device industries with the attitude of always making efforts.

HR MEDITECH C.E.O Hyeon-cheol, Lee





Dream of beautiful skin

HR MEDITECH is
a trueworthy company of renowned
plastic surgeons, and dermatologists in Korea.

HR MEDITECH PRODUCT


SONOCELLTriple	01
BODY MASTER	02
UV HIFU	03
POLKA RF	04
TRES	05
T-SONIC	06
COOL FIT	07
ENERGY UP	08
S-UP	09
ZERO PEEL	10
VITAL ION PLUS	11
MIUP HIFU	12
HIFU POINT	13
CELLUREON	14
MAGIC POWER	16
INDICELL	17
MAGIC CELL	18
MAGIC UP	19
MAGIC IMPACT · MAGIC POINT	20
BODY PLUS	21
DERMAWAVE	22
TRICELL · EXCEL LIPO	23
SKINPRO · SKINPROPLUS	24
MAGIC ICE · CRYOD	26
PLA-CLEAN	27
MAGIC PEEL · AQUA PLUS	28
MISTINANO · AQUA DUAL	29
ZIZIMI · AGUNG	30
CELLAXY AURORA	31
SUNLIGHT MASTER ⁶³³ · SELL LIGHT	32

Neat and clear in various modes

SONOCELL Triple

**Water droplet ultrasound,
firm and moist**


High-density ultrasound with wavelengths of 1, 3, 10, and 17MHz delivers energy to the entire skin layer.



[Various ultrasonic wavelengths
unique to SONOCELL Triple]

**Four wavelengths deliver sufficient
energy from shallow to deep layers
of the skin.**

With various wavelengths of 1, 3, 10 and 17MHz, SONOCELL Triple can care for the deeper layers of the skin, such as the muscle layer and the fat layer, and even care for skin troubles in the shallow layer.



[Handpiece composition and utilization]



• **3, 10, 17MHz (Blue LED)**

Sensitive skin care
Various skin trouble care
Helps calm skin irritation



• **1, 3MHz (Red LED)**

Affects the muscle layer and the fat layer
to relax the muscles, lymphatic massage,
skin swelling removal, care before lifting

A collection of endermologie, high frequencies, suction, and cavitation.

BODY MASTER

Perfect on the outside, line management

high frequency, suction, and color therapy LED Bipolar, ultrasonic wave, fat layer only. Selection and management cavitation, continuous with the addition of a contraction and relaxation system, Energy transfer, complete line care



[BODY MASTER One-Stop System]

Care for details regardless of wide or narrow parts!

Take care of it with a variety of functions!

BODY MASTER is cellulite (orange peel), bumpy fold, forearm, It takes care of detailed areas such as thighs and sides. It helps to make wannabe body shape with various functions.


[BODY MASTER effect]



▶ Before

▶ After

[BODY MASTER, Systematic line management]



ROLL-RF-VACUUM

APP-S-RF-VACUUM

LED BIPOLAR

CAVITAION

For your hidden v-line

UV HIFU

Perfect and confident elasticity

UV HIFU is managed from the inside without damaging the skin surface reborn with moisturizing skin texture as well as elasticity




[UV's high-intensity focused ultrasound]

Face and body management with high intensity focused ultrasound using the principle of a magnifying glass

Increased skin elasticity by applying high-intensity focused ultrasound to a specific area of the skin


[UV HIFU Function]

[UV HIFU effect]



▶ Before


▶ After



Deep folds of skin



UV HIFU Ultrasonic Energy Transfer



Without wrinkles

Fractional Bipolar

POLKA RF

Care for both epidermis and dermis

Use dedicated tips where desired
Epidermal Regeneration Effects
Transferring energy to the dermis.
Implementing the regeneration effect of high-frequency equipment




[POLKA RF Only Tools]

Heat stimulation to the skin layer to shrink aging collagen. Stimulates fibroblasts, helps induce new collagen production.


[POLKA RF Principle]

Perfect care from epidermis to dermis.
Enough energy transferable!

Fine dots, high-frequency energy.
Divided and minimal damage
Maximize management effectiveness



[POLKA RF Advantages]



The epidermal regeneration effect of Fractional and the high frequency dermis at once.
Helps induce collagen regeneration in the dermis.

Safe management with minimal management side effects

CAVITATION + RF VACUUM

TRES

Total management of cavitation, high frequency, and thermal suction with one device
 Simultaneously manage cellulite decomposition and skin elasticity
 Manageable in a short period of time and selectively manage the fat layer




[CAVITATION + RF VACUUM]

**Cellulite and Elasticity Management
 A combination for efficient total care!**

Manage by selecting only the fat layer due to ultrasonic waves A combination of cavitation and powerful suction high frequency.
 a sure before and after



[TRES Advantages]



Stimulate unnecessarily bloated fat by acting deep into the fat layer.
 Hard fat is softened and can be easily managed.

**Excellent effect on fat-rich
 abdomen, waist, hip, etc.**

A new paradigm for high-density ultrasound

T-SONIC

Elastic and moist skin with water droplets ultrasound


High-density ultrasound with wavelengths of 1, 3, 10, and 17Mhz delivers energy to the entire skin layer.



[Various ultrasonic wavelengths unique to T-SONIC]

Four wavelengths deliver sufficient energy from shallow to deep layers of the skin.

With various wavelengths of 1, 3, 10 and 17Mhz, T-SONIC can care for the deeper layers of the skin, such as the muscle layer and the fat layer, and even care for skin troubles in the shallow layer.



[Handpiece composition and utilization]



• **3, 10, 17MHz (Blue LED)**
Sensitive skin care
Various skin trouble care
Helps calm skin irritation



• **1, 3MHz (Red LED)**
Affects the muscle layer and the fat layer to relax the muscles, lymphatic massage, skin swelling removal, care before lifting

REDESIGN YOUR BODY

COOL FIT

COOL FIT's cooling system

For Cool Fit, after attaching the cool pad and handpiece to the area you want to be managed, it delivers cooling energy of minus 5-10 degrees to specific areas only.



Large Handpiece



Small Handpiece

[COOL FIT's cooling system]

**Stop the stressful diet anymore !
Start a comfortable diet with
COOL FIT's unique cooling system !**

COOL FIT delivers powerful cooling energy to the desired area such as the abdomen, sides, thighs, and arms to create the desired body shape.

[COOL FIT effect]




▶ Before



▶ After

[Changes in 90 days of using COOL FIT]



TO BALANCE YOUR BODY

ENERGY UP

Confirmed reliable care Rhythmic movement


In ENERGY UP, an energized manager receives a microcurrent using a living body current through a ground plate and delivers the microcurrent to a customer with a hand or handpiece.



[ENERGY UP Care Course]



The device transmits the microcurrent using the bio-current to the manager's feet through the ground plate.



The operator conducts a living body microcurrent to the ground plate.



The energized manager delivers the biological microcurrent to the customer by hand or handpiece.

[ENERGY UP effect]



▶ Before



▶ After

[Features of ENERGY UP]

Total care	Face and body total care possible by using biological microcurrent
Plates	Two managers can manage the upper and lower body with 4 ground plates
Handpiece	Two types of care systems that can be used with hand perfusion and pulse therapy systems.
Simple usage	Convenient management with simple usage

[PTS - 3 types of modes]

- MODE 1 : Care for local areas and lymph nodes
- MODE 2 : Systemic care for large areas, abdomen, etc.
- MODE 3 : Muscle care with multi-care mode

High frequency transmitted by fingertips

S-UP

Touch the skin

Unlike conventional high-frequency, multi-wavelengths are directly transmitted to the body through the S-UP globe.



[S-UP Grooming Care]

Clean and accurate management through S-UP gloves. Face & body all care possible !

Multi-wavelength reacts with the body and helps to convert cells in skin tissue into biothermal energy through movement of ions in the body and rotational movement of molecules.

[S-UP effect]



▶ Before

▶ After

[New concept multi-wavelength care system]



Body line care



Customized care



Comfortable care



Multi-wavelength



deep heat effect

PEEL + E.P + P.T.S

ZERO PEEL

Finding the answer to my own skin troubles.

Separates dead skin cells with 25KHz-28KHz ultrasonic vibration, and removes dead skin cells in two types, dry and wet.

Contour·dead skin cells·skin troubles simultaneously care

- * Ultrasonic vibration of 25KHz-28KHz per second
- * With a hook-shaped tip, it is possible to care for curved areas of the face
- * Ampoule penetration deep into the skin
- * Facial line care with muscle stimulation



[Handpiece composition and utilization]



①	PEEL	Removes dead skin cells, blackheads, and whiteheads, penetrates ampoule into the skin
②	P.T.S E.P	Use for face narrow area such as face
③	P.T.S E.P	Use for body wide area such as body

[ZERO PEEL effect]



▶ Before

▶ After

[PTS effect]

Prevent muscle atrophy / relieve muscle spasm
Facial lifting effect by stimulating the muscles of the face.
Effective on reflex nerve stimulation such as body/foot floor

Vitamin injector

VITAL ION PLUS

Penetrate high concentration of vitamin C into the skin

- 01 Use of ionized ampoule
- 02 Mode adjustable according to skin type
- 03 Skin dead skin cells, black head care
Nutrition and moisturizing care
- 04 Definite effect by administration of high concentration ampoule
- 05 Skin texture improvement and skin trouble relief



[VITAL ION PLUS principle]

Skin improvement by using the principle of pushing each other at the same pole and pulling the opposite poles together

Improving blemishes, various skin problems and aging skin by penetrating a large amount of high-concentration ionized vitamin C into the skin

[VITAL ION PLUS effect]



▶ Before

▶ After

[+, -, ± MODE usage tips]

*(-)After penetration of the ampoule due to the alkali reaction of the pole, it is acidified at the (+) pole.

*Vitamin ampoule administration after iontophoresis care before massage to remove oil from the skin

*Refrigerated storage of aqueous vitamin solution

[Lifting MODE usage tips]

*Care method that provides electrical micro stimulation to skin and muscles using a wet cotton swab

*Used for penetration of aging management ampoule

A new concept HIFU

MIUP HIFU

Reasons for flawless skin

MIUP HIFU is effective in skin elasticity and lifting care without damaging the skin surface.



[MIUP's high-intensity focused ultrasound]


Face and body management with high intensity focused ultrasound using the principle of a magnifying glass

Increased skin elasticity by applying high-intensity focused ultrasound to a specific area of the skin



*Ultrasound using the principle of a magnifying glass to treat only the desired area

[MIUP effect]



HIFU

POINT

Lifting without affecting everyday life

The HIFU system is a high intensity system that focused on ultrasound method of burning tissue using ultraviolet energy of 65-100°C, which occurs at the focus when ultrasound energy is collected in one place.


[Cartridge type]



[Product specification]

Type	D1.5	D3.0	D4.5
Frequency	7 MHz	4 MHz	4 MHz
Length	1.5mm	3.0mm	4.5mm
Life span	300,000 SHOT	300,000 SHOT	300,000 SHOT

[HIFU effect]



[Product standard]

Division	Standard
Rated Voltage	220V, 50/60Hz
Size	404.1 x 231.0 x 350.1(mm)
Weight	2.9kg

[Handpiece]




1. Operation indicator
2. Shot button
3. Cartridge

Heat Sensors

CELLUREON

New concept of decompression rolling system

Strong vibration generated by decompression capacity that is 15 times strong, makes rigid organ and muscles flexible. It is a new concept of decompression rolling system that facilitates blood circulation and promotes lymph circulation.



[CELLULION function]

Vacuum	a vacuum in the air inside the handpiece to reposition the fat and massage.
Suction	pulling caused by low pressure and vacuum. Technique and strength can be adjusted depends on degree of cup shape and suction.
Vibration	Generate heat inside the cell by repetitive movements periodically.
Acupressure	Make blood circulation and stimulate the nerve using nodules of the head.

[CELLULION body shaping]




It is a body management program that helps to manage cellulite by flexing rigid organ and muscles by improving circulation of blood and lymph by new concept of rolling system.

[CELLULION step 3]



[Thermotherapy]

Thermotherapy is a therapy that heals and improves the disease by maintaining some temperature of living organs or tissues in the range of 41°C - 45°C degrees within the limit to overcome the thermal effect while maintaining the physiological temperature control function.



► Before

► After



▶ 10.2 inch touchscreen



Flat type handpiece
- for use on back



Curved handpiece
- for use on arms and legs

Thermal
Vacuum



▶ Face handpiece




[Abdomen clinic once - 20min]



▶ Before



▶ After



▶ Before




▶ After

5 functions in one

MAGIC POWER

High frequency total care system

A system that combines the functions of RET, RF high frequency, Vacuum, LED, and Cryod into one.





[Handpiece name and function]



- ① VACUUM+RF: For Body
- ② VACUUM+RF: Small for body
- ③ VACUUM+RF: For Face
- ④ RET: For Body
- ⑤ RET: For Face
- ⑥ Scraping massage accessories
- ⑦ Scraping massage connector
*Separately purchased
- ⑧ RET-only tips
- ⑨ Cryod Handpiece

[LED RF high frequency Vacuum / RET high frequency]




RET high frequency
 Helps generate bioenergy deep in the human body.
 Helps regenerate elastin collagen fibersz
 Helps in fat dissolution and body line management

RET + BIPOLAR

INDICELL

Combine RET high frequency and bipolar!


RET is body shape management equipment that uses high-frequency bio-thermal energy in the body to generate deep heat from the dermis layer to accelerate body fat decomposition.



RET effect

Biologic energy is generated from the deepest.
 Rapid regeneration of rigid nodule masses and regeneration of elastin collagen fibers.
 It is effective not only for fat dissolution but also for body lifting.

[Handpiece name and function]



[Scraping massage]



▶ Scraping massage connector
 *Separately purchased

[INDICELL effect]



▶ Before ▶ After

MAGIC LINE+BIPOLAR

MAGIC CELL

Bipolar High Frequency + Vacuum + LED + Cryod

Efficient skin and body care with four functions of Bipolar High Frequency, Vacuum, LED, and Cryod in one device.



BIPOLAR
effect

- 01 Cell function activation
- 02 Body fat and body size reduction
- 03 Cellulite emission
- 04 Pain relief
- 05 Relieve muscle atrophy


[Product series]



MAGIC SLIM

MAGIC LINE

[LED Bipolar suction]



Characteristics

Four types of suction mode available
1 to 10 RF level adjustment function
Independently available function

[MAGIC ALL]



Vacuum

RF high
frequency

LED

Total vacuum system

MAGIC UP

Bipolar High Frequency + Thermal + LED + Vacuum


Bipolar High Frequency + Thermal + LED + Vacuum
Efficient skin and body care with 4 functions in one device



MAGIC UP Advantages

- * Bipolar high frequency + LED + vacuum, thermal + vacuum bipolar high frequency, vacuum, etc. can be used in various modes.
- * No plate
- * Independently available function
- * Integrated effect of RF high frequency and vacuum
- * Total vacuum system

[MAGIC UP effect]



Powerful vibration and decompression function and recuperative energy heat work simulataneously to give physical stimulation to the target area for muscle relaxation, blood circulation and lymph circulation promotion, cellulite, lipolysis, elasticity, facial contouring and lifting effects.

[3 in MAGIC UP]

Vacuum	The powerful 650 mmHg of physical force stimulates blood and lymphatic flow; promoting metabolism and facilitating the evacuation of wastes.
LED	Body light Therapy = 630nm LED Face light Therapy = 430nm LED
RF	When a multipolar high-frequency current of 1.0MHz is energized in the human body, molecules of the tissue vibrate each other as the direction of current changes and cause deep heat. Improving skin elasticity, promoting fat dissolution and emission

[MAGIC UP accessories]



► MAGIC UP accessories



► Vacuum accessories

Decompression rolling system

MAGIC IMPACT

Care of stiff muscles through the action of strong and arious vibration modes and heat, and help to promote blood circulation and lymph circulation.

[MAGIC IMPACT effect]

01 Vacuum **02** Pain relief **03** Muscle relaxation
05 Improve immunity **05** Fatigue recovery



► Before

► After / once 25min care

Energy heat + Thermal vacuum

MAGIC POINT

Various vacuum accessories can care for various parts of the body.

Various applications are possible due to the thermal vacuum function.



► Vacuum accessories




CET RF

BODY PLUS

Care system using high-frequency energy

Designed to effectively deliver high-frequency energy to the skin deep into the epidermis and sebum fat layer. It helps to improve the treatment area quickly.




[CET RF product series]



BARON


DIARO

[CET RF high-frequency]



It is a method of contacting the anode with the insulation coated ceramic and the grounded conductive plate. There is active electron activity centering on the ceramics contacted with the skin, and deep heat is generated. Therefore, it is effective in managing the dermis and subcutaneous fat close to the skin tissue.

[CET RF high-frequency effect]



- Lipolysis · normalization of cell function
- Hormonal balance reulation
- Promote lymph circulation
- Promotes contraction and regeneration
- Stimulate collagen layer
- Promote blood circulation

BIPOLAR

DERMAWAVE

Third-generation of high-frequency system

The bipolar high-frequency system is a third-generation of high-frequency system that has been stripped from the existing high-frequency ceramics coating system.


It is an advanced system that enables strong high frequency and efficient in skin care along with body management without coating.



BIPOLAR effect

- 01 Cell function activation
- 02 Body fat and body size reduction
- 03 Cellulite emission
- 04 Pain relief
- 05 Relieve muscle atrophy

[Product features]



- 01 There is no spark that cause fragile ceramics and skin damage.
- 02 Easy to operate, fast and powerful.
- 03 It stimulates the energy of fat cells by deep heat.
- 04 We offer two handpieces, so it is easy to concentrate on each site.

[Clinical date]



▶ FACE

▶ FACE & BODY

▶ BODY

▶ BODY

CAVI+BIPOLAR


TRICELL

Cavitation Frequency **58KHz**



▶ RF Handpiece

▶ Multi Handpiece




As the temperature and pressure in the bubble layer rises and the bubbles are destroyed, energy is generated, which selectively affects only the fat layer.

EXCEL LIPO

Find the answer to my own skin troubles

Ultrasound cavitation can be performed on all patients regardless of gender, and has excellent effects on obesity such as the fat-rich abdomen, waist, buttocks, and thighs.

[Egg burst simulation from cavitation power]



Comprehensive beauty machine

SKINPRO / SKINPROPLUS

All beauty care at once !

It is a multi-system with various skin care functions and various applicator for users can be provided for high treatment effect



- 01 Hygienic with single-use solution ampoule
- 02 Freely adjust solutions in aquafan
- 03 exfoliation, blackhead care, nutrition and moisturizing care at once
- 04 Reduced sebum, moist and clear skin tone
- 05 Simple and hygienic management

[SKINPRO]

Select only the desired applicator reconfigurable Minus Options System

*Request configuration of applicator without cryod - Possible



- 1. Tongs pole · 2. Grounding pole · 3. Iontophoresis
- 4. Scrubber · 5. Ultrasound · 6. Spray

[SKINPROPLUS]



- 1. Aqua peel
- 2. Cryod
- 3. Grounding pole
- 4. Tongs pole
- 5. Iontophoresis
- 6. Scrubber
- 7. Ultrasound

[Velvet mask penetration · Spray]



▶ Velvet mask penetration tong



▶ Spray

Comprehensive beauty machine



- 3 time modes
- Sebum exfoliation
- 4 stage energy level adjustment



- 3 time modes
- Regeneration due to vitamin infiltration
- Anti whitening facial cleansing
- 4 stage energy level adjustment



- Crease regeneration by collagen penetration
- Wrinkle lifting
- 4 stage energy level adjustment



- Provides 2 programs
- Pulse function control
- 4 stage energy level adjustment



- Velvet mask penetration tong
- It can be applied safely and effectively by giving physical stimulation to the target area.
- 4 stage energy level adjustment



[SMASTER]

[For export SKINPROPLUS]

[Product series]



5-in-1 Skin Care

- 01 Skin scrubber mode
- 02 Iontophoresis mode
- 03 Lifting mode
- 04 Velvet mask penetration mode
- 05 Ultrasonic mode



BARON



DIARO



Using the body's survival response

MAGIC ICE

CRYOD function and RF high frequency !

It is a principle to cleanse the body more cleanly by sending out unnecessary toxins and wastes in the body through the body's survival response.

MAGIC ICE
system

- 01 Powerful CRYOD therapy
- 02 Deep heat high frequency
- 03 Use of the two functions individually or together
- 04 Muscle relaxation
- 05 Cellulite reduction



- 1. Cold : -10°C / Warm : 55°C
- 2. Electroporation BODY
- 3. Electroporation FACE

CRYOD + HEATING + ELECTROPORATION

CRYOD

Electroporation

Electroporation applies an electrical pulse to the skin's surface to disrupt the skin's phospholipid membrane, temporarily creating a hole to penetrate the drugs needed to regenerate the skin.

Atopy
relief

Skin
soothing

Scalp
care

Pain
relief


Plasma

PLA-CLEAN

Plasma used in various places in life

Plasma has a sterilization effect and an infection prevention effect.

When plasma is generated, harmful molecules around it are decomposed to remove harmfulness.




[Plasma]

Plasma is ionized by ion separation into electrons and atomic nuclei.

Plasma exhibits excellent effects by effectively destroying cell membranes due to simultaneous physical and chemical action.

Staphylococcus aureus, which is the main cause of atopy and itching, has a strong resistance and is not easily eradicated, but experiments that have been killed after plasma treatment have been published in various papers.

[Plasma Staphylococcus aureus killing test]



[PLA-CLEAN effect]



▶ Before

▶ After / once 30min care



▶ Before

▶ After / once 30min care



Supply plenty of moisture and
remove the waste of the face

MAGIC PEEL

Moisturize, remove face wastes!

AQUA PLUS PEEL uses a three step of special solution to inhale acne, blackhead, whitehead, exfoliation, folliculitis, dust, fine dust, etc. as the solution injected from the center of headpiece tip rapidly circulates based on the scientifically designed handpiece tip.

MAGIC
PEEL

- 01 Remove sebum, acne
- 02 Moisture supply
- 03 Improvement Skin texture
- 04 Skin peeling·Folliculitis reduction



Electroporation



Skin care through water particles

AQUA PLUS

A system used for skin cleansing and cosmetic application by spraying moisture particles on the skin and scalp with powerful pressure

AQUA
PLUS

- 01 Nutrient penetration into skin
- 02 Available all facial scalp treatment
- 03 Spray jet intensity control
- 04 Pressure adjustable

LOW LEVEL LASER + LONTOPHORESIS
NANONIZATION MIST + EMS

MISTINANO

Multi-manage with MISTINANO

MISTINANO can manage the whole body with four functions! beauty device that can work simultaneously with lifting + whitening + anti-inflammatory + toxin emission + EMS.



► Before

► After



Sanitary and powerful

AQUA DUAL

AQUAPEEL manage easier and simpler

It's easy, hygienic, and easy to adjust the compression directly from the pen. Powerful spray mist to help with ampoule penetration.



► Before

► After



AQUA DUAL
solution

Far Infrared radiation steamer

ZIZIMI / AGUNG

If your body temperature goes up 1 degree, your immune system will increase 5 times.


Far infrared

It refers to infrared rays with a wavelength of 25mm or more. it is invisible because it has a longer wavelength than visible light. Also has a large heat effect and has strong penetration. Among these infrared rays, the longest wavelength is far-infrared.

[source : Doosan encyclopedia]



[Effect on heat action]



▶ Before


▶ After

it is effective to prevent various diseases such as aging, metabolism promotion, and chronic fatigue by activating the cell tissue by shaking the cell minutely 2,000 times a minute. When it comes into contact with moisture and protein molecules that make up the cell.

[Excellent Far-infrared emissivity]

[ZIZIMI Product standards]

Division	Standards	
Rated Voltage	220V, 60Hz	
Size	FLAT(mm)	510 × 80 × 360 (H)
	BEND(mm)	510 × 490 × 240 (H)
Weight	4.0 kg	



Excellent Far-infrared emissivity of 89.6% at 90°C


The Far-infrared emissivity rises in proportion to the temperature, therefore the emissivity at the actual operating temperature is important.

7-in-1 PDT Ca+ Light

CELLAXY AURORA

Seven Happinesses as One


The PDT color therapy has different effects depending on the frequency of each color. Adjust the calcium ion concentration to suit your own skin concerns.



[Product features]

- 01** Easy calcium filling and color therapy
- 02** GUI with high readability and clean design
- 03** Touch lighting in the dark OK.
- 04** Intensive care for areas that are easy to miss, such as forehead, under eye, mouth, etc.

[Frequency by Color]



[Color-specific effects]

RED 635nm	Pain relief, blood circulation promotion skin elasticity enhancement
GREEN 560nm	Induce sensitivity, pain relief collagen production
BLUE 415nm	Pore contraction, acne sterilization skin trouble soothing
PURPLE 308nm	Lymphatic management, waste discharge, skin tissue activation
YELLOW 590nm	Skin whitening, skin tone-up immune enhancement
ORANGE 610nm	Metabolic, blood circulation promotion, sensitive skin




Whitening Management Program


SUNLIGHT MASTER⁶³³

Use German original 633nm wavelength lamp

- 01 Bluetooth function eliminates boredom during management
- 02 Promote collagen + elastin synthesis
- 03 Real 633nm wavelength lamp in Germany penetrates deep into the skin
- 04 Effect on dull skin, aging skin, acne, and troubled skin



► Before



► After




- 1. Tornado vibration massage
- 2. Cold : -10°C / Warm : 55°C)

Tornado vibration massage

SELL LIGHT

The tornado vibration massage system which causes 8000 waves per minute

Is effective in improving wrinkles and skin elasticity.



MEDICAL DEVICE

OEM · ODM

BEAUTY DEVICE

SKIN CARE DEVICE



17, Janggogae-ro 231beonan-gil, Seo-gu, Incheon, Republic of Korea
T. +82. 32. 713. 8008 | F. +82. 2. 532. 4673 | H. www.hrmeditech.co.kr